

STARLAB focus

STARLAB Customer Information



A technical breakthrough!



TipOne RPT

STARLAB's new *Repel Polymer Technology (RPT)* optimises conventional polypropylene surfaces to improve the flow of "tricky" liquids like various biological buffers, Triton-X or glycerin. STARLAB's TipOne® RPT have no coating to interfere with your sample and offer an extremely hydrophobic surface.

In contrast to other technologies that use siliconised or fluoro-coated surfaces, TipOne® RPT tips offer exceptional delivery without the risk of interaction with additives or treatments, even after sterilisation or autoclaving.










Features

- No coating! Made from virgin polypropylene
- Improved sample recovery
- Reduced DNA binding
- Autoclavable
- Fits all popular pipettes
- Certified RNase, DNase, DNA & Pyrogen free



TipOne RPT



ORDERING INFORMATION	RACKS	STERILE RACKS	REFILLS
	10 x 96-tip racks	10 x 96-tip racks	10 x 96 refills
 10 µl RPT Graduated Tip (graduations: 2.5 µl, 10 µl)	n/a	n/a	S1161-3700
 10 µl XL RPT Tip (graduation: 10 µl)	S1160-3800	n/a	S1160-3700
 200 µl RPT Bevelled Tip	S1161-1800	n/a	S1161-1700
 1250 µl XL RPT Tip	S1161-1820	n/a	S1161-1720
 10 µl XL RPT Filter Tip (Sterile) (graduation: 10 µl)	n/a	S1180-3810	n/a
 20 µl RPT Bevelled Filter Tip (Sterile)	n/a	S1180-1810	n/a
 100 µl RPT Bevelled Filter Tip (Sterile)	n/a	S1180-1840	n/a
 200 µl RPT Graduated Filter Tip (Sterile) (graduations: 50 µl, 100 µl, 200 µl)	n/a	S1180-8810	n/a
 1000 µl XL RPT Filter Tip (Sterile)	n/a	S1182-1830	n/a

All tips shown actual size. n/a = not available

STARLAB GmbH

Kurt-Fischer-Straße 22
D-22926 Ahrensburg
Germany

T: +49 (0)4102 60440
F: +49 (0)4102 604420
info@starlab.de

STARLAB (UK), Ltd

4 Tanners Drive
Milton Keynes MK14 5NA
United Kingdom

T: +44 (0)1908 283800
F: +44 (0)1908 283802
infoline@starlab.co.uk

STARLAB-FRANCE

38 Avenue Henri Barbusse
92220 Bagneux
France

T: +33 (0)1 45 3652 80
F: +33 (0)1 45 3604 86
info@starlab-france.com

STARLAB Srl

Via Corelli 36/5
20134 Milano
Italy

T: +39 02 7020 1040
F: +39 02 7020 1033
info@starlabitalia.it

STARLAB Scientific India Private Ltd.

First Floor, Old No: 322, New No: 24
Thambu Chetty Street
Chennai - 600001, India

T: +91 (0) 44 42130702 / 703
F: +91 (0) 44 42130701
info@starlab.co.in

www.starlabgroup.com

